

...for multiple functions



GRADUATED COMPRESSION

- Supports articulations where necessary
- produces an intermittent **massage** which, in combination with the movement of the muscles, stimulates blood and lymph circulation
- the stimulated **circulation of the tissues** promotes their oxygenation and nourishment, accelerating the healing process
- stimulated **lymphatic circulation** prevents the formation of oedemas or promotes their reabsorption



PROPRIOCEPTIVE EFFECT

facilitates control of movements, reducing pain



MAXIMUM WEARABILITY, COMFORT, AND FREEDOM OF MOVEMENT

9000P00199INT01.0324

lastik PLUS

M90132F - Epilastik Plus Elasticated fabric elbow support
M90119FR - M90119FL - Manulastik Plus Elastic wrist support
M90104F - Genulastik Plus Mesh elastic knee brace
M90137F - Malleolastik Plus Elastic knit ankle support

Code	EPILASTIK PLUS M90132F						
Size	0	1	2	3	4	5	6
Forearm circumf. 13 cm from the epicondyle cm	17/19	19/21	21/23	23/25	25/27	27/29	29/31
Color	yellow fluo						

fits right and left

Code	MANULASTIK PLUS M90119FR - M90119FL						
Size	0	1	2	3	4	5	6
Wrist circumf. cm	13/14	14/15	15/16	16/17	17/18	18/19	19/20
Color	yellow fluo						

specify right or left

Code	GENULASTIK PLUS M90104F							
Size	0	1	2	3	4	5	6	7
Thigh circumf. 14 cm above the middle of the kneecap cm	35/38	38/41	41/44	44/47	47/50	50/53	53/56	56/59
Calf circumf. 12 cm below the middle of the kneecap cm	25/28	28/31	31/34	34/37	37/40	40/43	43/46	46/49
Color	yellow fluo							

fits right and left

Code	MALLEOLASTIK PLUS M90137F						
Size	0	1	2	3	4	5	6
Ankle circumfer. cm	15/17	17/19	19/21	21/23	23/25	25/27	27/29
Color	yellow fluo						

fits right and left



Indications


Epilastik Plus

- Lateral epicondylitis (tennis elbow) and medial epicondylitis (golfer's elbow)
- Chronic pain from bursitis or tendinitis



Manulastik Plus

- Wrist and carpal sprains
- Tendinitis wrist · Arthrosis and arthritis of wrist
- After wrist fractures
- Post-surgery wrist treatment



Genulastik Plus

- Knee pain due to bruises
- Early stage arthrosis
- Chondropathy/ Chondromalacia patella
- Mild patella instability
- Prevention of patella luxation during sports.



Malleolastik Plus

- 1st degree tibiotarsal sprains
- Peroneal or instep tendinopathy

lastik PLUS

Love sport,
take care
of your body



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



orthoservice.com



ORTHOSERVICE
RO+TEN
Take care feel better



ORTHOSERVICE
RO+TEN
Take care feel better

The best of features...



DYNAMIC MESH FABRIC

having 4 different weaves with different holds:

- **MEDIUMPRESSION** provides "medium" compression in the supramuscular areas
- **STRONGPRESSION** provides firm compression in the peri-articular areas
- **WRINKLE RAIL** accompanies the flexion-extension of the joint
- **SOFTPRESSION** that avoids a "tourniquet" effect in the areas that come in contact with the skin

... for unique performances



MICRO SILICON TUBE

for more **softness** on the skin and a **reinforcement** of the orthosis while its being worn



ANATOMICAL DESIGN 100%

to ensure **effectiveness** and maximum comfort



epilastilk PLUS

manulastilk PLUS

genulastilk PLUS

malleolastilk PLUS



SILICONE INSERTS

to **relieve** compression around bony prominences